## GEN 4: YIELD CALIBRATION





3



Select 'Harvest', then 'Yield Calibration'.

5

Select 'Start'.



Select 'Machine Settings', then 'Calibrations & Procedures'.



## Select 'Record'.

2

4

6



Once you have collected a sufficient sample, select 'Done'.



## GEN 4: YIELD CALIBRATION

8



9

HORIZON

## Select 'Match Scale Weight'.



Enter the scale weight of the sample load, then select 'Save'.

Calibrations & Procedures 🔞								
Sample   Enter Scale Weight								
Recorded Weight	20007							
30,424 🌲								
Sample Details	7 8 9							
Date/Time Collected: Flow Rate Category: Mass Flow Rate Avg: ~0.0 lbs/sec	4 5 6							
Duration: 0 minutes Field:								
	+/- 0 .							
📋 Delete Sample 🗶 Cancel 🗸 Save								

Select the sample you'd like to match.

🛱 Calibrations & Procedures 🧃						
Match Scale Weight 🥡						
Sample ID	Quality	Date / Time	Sample	Scale	Margin	
0001			30,424 lbs			
0002			35,900 lbs			
0003			15,820 lbs			
		1-				